



Knee Support



Owners Manual and Instructions for Use



**IMPORTANT
READ THIS
MANUAL
BEFORE USE**

Warranty

This product is guaranteed against all defects in manufacture for a period of one year.

This does not affect your statutory rights

Warranty number:

Operating Instructions

Knee Support

Intended use:-

The Knee Support is intended as a comfort aid for the client, it is not intended to provide postural support.

Clients should be able to support themselves independently in a sitting position in the bath.

Using your Knee Support

Ensure the bath is clean and there is no oily or soapy residue on the bath surface.

The Knee Support is fitted to the bath surface using suction cups. Choose the position where the knee Support is to be fitted then place the Knee Support in the bath and press firmly down to ensure it is fitted securely.

We advise against the use of oil based products when using a knee Support.

Cleaning

The Knee Support may be cleaned using most mild cleaning agents. The Knee Support should be cleaned using a warm soapy solution and rinsed with a damp cloth. Attention should always be paid to instruction labels on the outside of cleaning agent containers.

Chemical agents must not be used, i.e. drain cleaners, acetone, oven cleaners or acids.



**KNEE SUPPORT ONLY.
MUST NOT BE USED
AS A BODY BRACE**

Troubleshooting

In the unlikely event that problems occur with the unit, carry out the following easy

Fault: The Knee Support does adhere to the bath surface

Cause: There is a soap or oily residue on the bath surface

Remedy: Clean the bath surface using a mild cleaning agent and rinse thoroughly with clean warm water.

If the problem persists after carrying out these remedies, please contact Astor-Bannerman using the details on the back cover.



Unit 11f, Coln Park
Andoversford
Cheltenham
Glos, GL54 4HJ
Tel: 01242-820820
Fax: 01242-821110

e mail: sales@astorbannerman.co.uk